

# Lansetter och blodprovstagare

## Unistik® 3 – Engångslansetter med Comfort Zone®-teknologi.



**Gauge:**  
28 G  
**Depth:**  
1,8 mm  
**Blood Volume Range:**  
5 µL – 30 µL



**Gauge:**  
23 G  
**Depth:**  
1,8 mm  
**Blood Volume Range:**  
30 µL – 75 µL



**Gauge:**  
21 G  
**Depth:**  
2,0 mm  
**Blood Volume Range:**  
75 µL – 125 µL



**Gauge:**  
18 G  
**Depth:**  
1,8 mm  
**Blood Volume Range:**  
200 µL +

Varunummer	Produktbeskrivning	Antal/förp.
203825	Unistik 3 – Comfort	100
AT1044	Unistik 3 – Comfort	200
205413	Unistik 3 – Normal	100
AT1004	Unistik 3 – Normal	200
205412	Unistik 3 – Extra	100
AT1052	Unistik 3 – Neonatal	100

## Blodprovstagare för egenvård - Autolet® Impression

Kan förskrivas på hjälpmedelskort och passar alla lansetter till stickhandtag (utom ACCU-CHEK® Softclix). Autolet® har Comfort Zone®.



Varunummer	Produktbeskrivning	Antal/förp.
202483	Autolet® Impression	1

SE-1020-140822