

Lansetter

Unistik® 3 – Engångslansetter med Comfort Zone®-teknologi.



Gauge:
28 G
Depth:
1,8 mm
**Blood Volume
Range:**
5 µL – 30 µL



Gauge:
23 G
Depth:
1,8 mm
**Blood Volume
Range:**
30 µL – 75 µL



Gauge:
21 G
Depth:
2,0 mm
**Blood Volume
Range:**
75 µL – 125 µL



Gauge:
18 G
Depth:
1,8 mm
**Blood Volume
Range:**
200 µL +

Varunummer	Produktbeskrivning	Antal/förp.
203825	Unistik 3 – Comfort	100
AT1044	Unistik 3 – Comfort	200
205413	Unistik 3 – Normal	100
205412	Unistik 3 – Extra	100
AT1052	Unistik 3 – Neonatal	100